



Fresenius Medical Care North America Service Bulletin


Equipment: All 2008 Series Hemodialysis Machines
Bulletin: 18-FRS-001 Rev A
Subject: New UF Pump Filter Part Number & Appearance

1.0 PURPOSE


To inform all Field Service and Technical Support personnel as well as Fresenius Medical Care-trained customers of the new part number and appearance of the UF Pump Filter used in all 2008 Series hemodialysis machines:

2.0 SCOPE


The UF Pump Filter (74) has a new part number (P/N **M48341**) and appearance. Even though the appearance has changed, the specifications for the new filter remain the same as the old filter; therefore, these UF Pump Filters are interchangeable.



Old UF Pump Filter
(P/N 672574)



New UF Pump Filter
(P/N **M48341**)



UF Pump Filter (74) Hydraulic Location

If there are any questions regarding this bulletin, contact Fresenius Medical Care Technical Support at 800-227-2572.

© 2018, Fresenius Medical Care North America, All rights reserved.