



**FRESENIUS  
MEDICAL CARE**

## Fresenius Medical Care North America Service Bulletin


Equipment: GranuFlo® Dissolution Tank - 132 Gallon  
Bulletin: 13-GRN-001 Rev A  
Subject: Tank Twist Lid and Ring

### 1.0 PURPOSE

To inform all Field Service and Technical Support personnel of a new Tank Twist Lid and Ring Assembly for the 132 gallon GranuFlo® Dissolution Tank.

### 2.0 SCOPE

A new assembly (P/N 337002-12) has been created for the 132 gallon GranuFlo® Dissolution Tank. It combines the tank twist lid and the twist lid ring into one assembly. The figure below shows both parts that make up the new assembly:



Tank Twist Lid and Ring Assembly (P/N 337002-12)

The tank twist lid (P/N G84-001-12) and the tank twist lid ring (P/N G84-001-13) are no longer available as separate spare parts.

If there are any questions regarding this bulletin, contact Fresenius Medical Care Technical Support at 800-227-2572.

© 2013, Fresenius Medical Care North America, All rights reserved.