

# 2008T BlueStar Hemodialysis Machine In-Service Training Checklist



Facility: \_\_\_\_\_ Date of In-Service Training: \_\_\_\_\_

Address: \_\_\_\_\_ Contact: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_ Phone: \_\_\_\_\_

Presented by: \_\_\_\_\_, RN

## Training Topics

Topics may be modified to accommodate the specific policies and procedures of the unit being trained. Check marks indicate topics reviewed during the In-Service training; no representation is made regarding the participant competence or expertise.

Training Topics	✓
<b>1. 2008T BlueStar Features</b>	
• Familiarize Physical Features	
• Review of Screen Attributes	
• Review of Alarms/Troubleshooting	
• Disinfection & Maintenance	
• Auto Start/Testing	
• Sustained Low Efficiency Dialysis (SLED)	
• Low Volume Mode	
• Disinfect Log	
• Auto Prime	
• Assisted Reinfusion	
• Independent Conductivity/pH Testing	
• PatientCard System	
– Entering Patient Information	
– Default Parameters	
– Dialysate/Heparin/BP/UF Settings	
– Prescription	
• Alarm Volume Controls	
• Select Program Screen Changes	
- DIASAFE® Plus Reminder	
- PM Reminder	
- Digital Hour Meter	
- Digital Serial Number	
• New Heparin Syringe Choices	
• Conductivity Alarm Position Lock (if clinic is using)	
• New Treatment Prescription Warnings	

	✓
• Low Power Mode	
• Idle Mode Option of 100 ml/min. Dialysate Rate	
• New 150 ml/min. Dialysate Rate	
• Isolated UF	
• Alarm Volume Control, includes: Difference Between Low and High Priority Sounds, Acute vs. Standard (Chronic), How to Adjust Volume, Mute Symbol, and Volume Control Limits.	
<b>2. Service Mode Options</b>	
• Default Dialysate Flow 500/600/700/800 ml/min.	
• Plasma Sodium — Hide/Show	
• Sodium Variation System (SVS) On — Yes/No	
• Idle Mode Dialysate Flow Rate — 100/300 ml/min.	
• Conductivity Alarm Position — Lock/Unlock	
• New Treatment Prescription Warning	
• Default Rx Screen	
• Low Volume Max UF Rate — 500 – 1000 ml.	
• Maximum UF Per Treatment vs. Per Hour	
• Alarm Type — Acute/Chronic	
• Maintenance Scheduler Options	

# 2008T BlueStar Hemodialysis In-Service Sign In



Print Name	Signature	Title

